



Assistive Device Catalogue

Svastya – a health initiative by GiftAble!

Crutch

Crutch is a mobility aid that transfers the body weight from lower body (the legs) to the upper body (hands & shoulders). It is useful for people who cannot use their leg to support their body weight due to short-term injuries, neurological impairment or disability.



Note: The color and design of these product may vary from what is shown in the images.

Crutch



Single Crutch

Features:

- Under arm padded
- Aluminum frame
- Adjustable hand grip
- Adjustable length
- Rubber pod



Double Crutch

Features:

- Under arm padded
- Aluminum frame
- Adjustable hand grip
- Adjustable length
- Rubber pod

Crutch



Elbow Crutch

Features:

- Anatomical Palm Grip
- Aluminum frame
- Adjustable length
- Light Weight
- Rubber Pod



Tripod/Quadripod Crutch/ Stick

Features:

- Anatomical Palm Grip
- Aluminum frame
- Adjustable length
- Light Weight
- Rubber Pod

Walker

Ambulation is the ability to walk from place to place, but this can be difficult for children and adults with physical disabilities. A walker enables them to be mobile and perform their daily activities to their maximum functional capacity.



Note: The color and design of these product may vary from what is shown in the images.

Walker

Adjustable Walker



Features:
Platform Walker,
Height Adjustable
Aluminum frame
Seat width 47 cm
Light Weight
Rubber Pod

Foldable Walker



Features:
Platform Walker,
Height Adjustable
Aluminum frame
Seat width 47 cm
Light Weight
Easy to carry
Rubber Pod

Walker

Cerebral Palsy Walker



Features:
wheels that are 8" in diameter.
Larger wheels & handle
Rough terrain walker
Light weight

Rollator



Features:
Platform Walker,
Easy to move
Light Weight
Rubber Pod/ wheel

Orthosis

Orthoses are custom-made or prefabricated devices used to optimize body functions and body posture to help promote an individual's participation in everyday activities and therefore the community. Using this assistive device also prevents secondary complications such as deformities.



Note: The color and design of these product may vary from what is shown in the images.

Orthosis

Ankle Foot Orthosis



Features:
Light weight
Easy to wear
Can fix within shoe

Knee Ankle Foot Orthosis



Features:
Light weight
Easy to wear
Can fix within shoe
Knee lock
Well padded
Available with shoe & without shoe

Orthosis

Supra Malleolar Orthosis



Features:
Light weight
Easy to wear
Can fix within shoe
Velcro strapping
Well padded

Leaf Spring Orthosis



Features:
Light weight
Easy to wear
Can fix within shoe
Well padded
Can be customize

Splint

A **splint** is defined as "a rigid or flexible device that maintains in position a displaced or movable part, also used to keep in place and protect an injured part" or as "a rigid or flexible material used to protect, immobilize, or restrict motion in a part".

Splints can also be used to relieve pain in damaged joints.



Note: The color and design of these product may vary from what is shown in the images.

Splint

Aeroplane Splint



Features:

Light weight
fixed by long M. S.
strip with Horizontal
Palm support Rubber
band & leather loop
for thumb
Well padded

Cock Up Splint



Features:

Light weight
fixed by long M. S.
strip with Horizontal
Palm support Rubber
band & leather loop
for thumb
Well padded

Chair/ Frames/ Board

A frame is an assistive device that is specifically designed to allow individuals to stand or sit from short to extended periods. Using a frame aids the development of motor coordination and head control. Also facilitates muscle and bone development. This increases the individuals' self-confidence and allows them to experience the joy of standing and sitting.



Note: The color and design of these product may vary from what is shown in the images.

Chair/ Frames/ Board

Cerebral Palsy Chair



Features:
Wooden frame
Cushion
Pudum
Seat Adjustable

Corner Chair



Features:
Customize design
Easy to clean
Chest and head strap
Well padded
Pudum

Chair/ Frames/ Board

Standing Frame



Features:

- Customize design
- Easy to move
- Chest and head strap
- Well padded
- Pudum

Scooter Board



Features:

- Hard plastic
- 360 degree can rotate
- Front and back movt.

Gaiters/ Foot wear

Gaiters



Features:
Made of cloth
Customize design
Easy to wear

CTEV Shoe



Features:
Customize design
Well padded

Prosthesis

Below Elbow Prosthesis



Features:
Battery operated
Customize design
Easy to wear

Above Elbow Prosthesis



Features:
Battery operated
Customize design
Easy to wear

Prosthesis

Below Knee Extension Prosthesis



Features:
Customize design
Easy to wear

Above Knee Extension Prosthesis



Features:
Customize design
Easy to wear

Prosthesis

Symes Prosthesis



Features:
Customize design
Easy to wear

Brace

Ash Brace



Features:

- Aluminium body.
- Etha foam covered pads.
- Better impact absorption
- Enhanced comfort,
- Reduced pressure on pubis and chest.
- Anatomical and rigid back pad.
- Customizable.

Lumbar Sacral Brace



Features:

- Flexible sizing.
- Double pull mechanism.
- Well cushioned back pad.

Brace

Taylor Brace



Features:

- Easy to wear
- Customize as per age
- Provides free shoulder movement
- No obstruction of the thoracic

Boston Brace



Features:

- Customize design
- Easy to wear
- prefabricated pads

Collar

Cervical Collar



Features:

- Easy to wear
- Latex free foam
- Size can be adjustable

Philadelphia Collar



Features:

- Non-toxic hypo allergenic
- Latex free foam
- Easy to wear
- Size can be adjustable

Hearing Aids

A hearing aid is a device designed to improve hearing by making sound audible to a person with hearing loss.



Hearing Aids

A&M XTM P A4



Features:

- 8 channel
- Fit to Go App
- Rocker switch
- Compactable with smart remote app

Signia Run SP BTE



Features:

- 8 channel
- Feedback Cancellation
- Rocker switch
- Signal to Noise ratio 24dB SPL

- Cost shall vary depend upon model, brand and condition
- Moulding and maintenance service provided
- Speech therapy services available

Hearing Aids Brands



Wheelchair

Motivation has been at the forefront of wheelchair design and provision. All these wheelchairs are fully adjustable so they can be properly fitted to the user, reducing the risk of life-threatening secondary complications such as pressure ulcers.



Note: The color and design of these product may vary from what is shown in the images.

Wheelchair

Moti-Start



Features:

From 6 month to 2 yrs age group

ADL board

Adjustable seat depth, backrest height

Size can be adjustable

Wheelchair

Moti-Go



Features:

From 3yrs to 17yrs age group

Three-wheels to provide stability and easy riding over uneven ground

ADL board

Adjustable seat depth, backrest height

Available in both medium and large size

Wheelchair



Children Activity Wheelchair

Features:

- Three-wheels to provide stability and easy riding over uneven ground
- Adjustable seat depth, backrest height
- Removable rear wheels with a durable quick release mechanism for easy transportation

Other services



Fostering Inclusion



Therapeutic Intervention

- Physiotherapy
- Occupational Therapy
- Special Education
- Early Intervention
- Speech, Play & Group Therapy
- Music Therapy
- Art Based Therapy
- Yoga
- Counseling



School Program

- Sensitization of National Education Policy (NEP)
- Globally Accredited
- Standardize Assessment
- Teacher training
- Parent training
- Peer training
- Therapy



Social Security Scheme

- UDID card
- NIRAMAYA- Medical Insurance
- State Govt. Pension linkage



**For more information about our products,
or to find out about training package to suit your needs,
please contact us!**

Email ID: nirakar.devdas@giftabled.org / divya.s@giftabled.org

Contact No: +91 9880494953 / +91 9776637201

Website: www.svastya.giftabled.org

Address: #148, 1st A main East of NGEF layout, Kasturi Nagar, Bangalore 43.